



QUICOLOR NF[®]

RAPID ANTIBACTERIAL SUSCEPTIBILITY TESTING MEDIUM
FOR NON-FERMENTATIVE BACTERIA



Catalogue #: MP 0250

Instructions for Use

For In Vitro Use Only

Product's name:

QUICOLOR NF[®]

Product's intended use:

QUICOLOR NF[®] is a rapid antibacterial susceptibility-testing medium for non-fermentative bacteria like *Pseudomonadaceae* and *Acinetobacter* species. It is intended for **in vitro use only**.

General information:

Methods of disk diffusion and dilution in liquid or solid media are widely used to determine the susceptibility of bacteria to antibacterials. These methods are time consuming and require overnight incubation.¹

Early determination of antibacterial susceptibility of bacteria, isolated from patients with meningitis, bacteremia and sepsis, is very important for the selection of convenient therapy as soon as possible. Even when there is no such urgent need, determination of antibacterial susceptibility on the same day the culture results are obtained, has many advantages such as choosing the appropriate treatment regimen early, lowering the number of side effects, decreasing the duration of stays in hospital and the amount of expensive antibacterials used and thus total health care costs.²⁻⁵ The color of **QUICOLOR NF[®]** which is yellow, changes to red by the metabolic activity of bacteria growing on the agar surface. Inhibition zones appear in 4 to 6 hours, much before the bacterial layer forms.

Limitations of the method:

Some slow growing species of non-fermentative bacteria may infrequently require more than 6 hours for the completion of the test.

Principles of the procedure:

In classical Kirby Bauer disk diffusion antibacterial susceptibility testing, the inhibition zones are visualized when a bacterial layer on the surface of agar forms usually at 18 to 24 hours⁵. **QUICOLOR NF[®]** changes its color due to metabolic products of bacteria growing on the surface, which occurs within 4 to 6 hours, much before the bacterial layer forms. The inhibition zones around the disks containing antibacterials stay at the original color of the media. Thus

inhibition zones appear as yellow zones in red areas where bacteria grow.

Ingredients:

Peptones, carbohydrates, salts, vitamins, color dye indicators.

Cautions and warnings:

- Always use freshly cultured bacteria to obtain susceptibility results within 4 to 6 hours.
- Be sure that the tested strain is obtained as pure culture.

General safety precautions:

- FOR *IN VITRO* USE ONLY.
- Always wear gloves when working with potential biohazard material.
- Never use mouth pipetting.
- If spills of the contaminated material occur, disinfect with 2.5% hypo-chlorite solution.
- All contaminated material should be discarded according to biosafety rules.
- At a minimum, universal precautions should be employed.

Storage instructions:

Store at 2-8°C, protected from direct light.

Indications of instability or deterioration:

QUICOLOR NF[®] should not be used if there are signs of contamination or deterioration (discoloration, drying, shrinking, cracking).

List of materials provided:

MP 0250: Ten **QUICOLOR NF[®]** 90mm agar plates in each pack.

MP 0250-12: Ten **QUICOLOR NF[®]** 120mm agar plates in each pack.

MP 0250-15: Ten **QUICOLOR NF[®]** 150mm agar plates in each pack.

Additional material required:

- Sterile suspension fluid.
- Wireloop.
- Bunsen burner.
- Sterile swab.
- Antibiotic disks.



QUICOLOR NF®

Instructions for use:

1. Prior to inoculation, the media should be brought to room temperature.
2. Prepare a bacterial suspension of 0.5 to 1.0 Mc Farland from a fresh culture.
3. Spread the suspension onto **QUICOLOR NF®** using a sterile swab.
4. Place the paper disks containing the antibacterials to be tested on the surface of **QUICOLOR NF®**.
5. Incubate at 35-37°C.
6. The colored inhibition zones appear within 4 to 6 hours.
7. Susceptibility is determined by measuring the diameter of the inhibition zone.

Evaluation of the results:

The evaluation of the results should be done according 'QUICOLOR® disk diffusion, inhibition zone interpretive standards' provided together with this package insert. These standards are obtained by comparison of the susceptibility results of a large number of bacterial isolates obtained from clinical samples and studied by Kirby-Bauer disk diffusion method.⁵

Quality control:

Positive control: *Pseudomonas aeruginosa* ATCC 27853.

Negative control: Uninoculated plates.

Time restrictions:

The colored inhibition zones should be evaluated when they form at 4 to 6 hours. If the zones cannot be evaluated within 8 hours then the color zones will disappear. However inhibition zones will reappear after overnight incubation in the form of bacterial layer on the surface as in classical disk diffusion method.

The following may impede the visualization of inhibition zones within 4 to 6 hours:

- 1- Preparation of the inoculum from an old culture plate. (The inoculum should be prepared from one day old culture that is not kept in refrigerator.)
- 2- Preparation of the inoculum with turbidity less than 0.5 Mc Farland.
- 3- Leaving uninoculated parts on the medium when spreading the inoculum by the swab.

Limitations of the procedure:

Antibacterial susceptibility results may not be obtained within 4 to 6 hours if bacteria to be tested are not obtained from fresh culture plates.

Performance characteristics:

When susceptibility test using **QUICOLOR** is performed according to instructions, over 95% overall agreement should be expected between the results of Kirby Bauer and **QUICOLOR** disk diffusion methods.⁶

Bibliography:

- 1- Barenfanger, J., C. Drake, G. Kacich. 1999. Clinical and financial benefits of rapid bacterial identification and antimicrobial susceptibility testing. J. Clin. Microbiol. 37:1415-1418.
- 2- Doern, G., R. Vautour, M. Gaudet, and B. Levy. 1994. Clinical impact of rapid in vitro susceptibility testing and bacterial identification. J. Clin. Microbiol. 32:1757-1762
- 3- Schifman, R., A. Pindur, and J. A. Bryan. 1997. Laboratory practices for reporting bacterial susceptibility tests that affect antibiotic therapy. Arch. Pathol. Lab. Med. 121:1168-1170
- 4- Granato, P. A. 1993. Impact of same day tests versus traditional overnight testing. Diagn. Microbiol. Infect. Dis. 16:237-243
- 5- Performance Standards for Antimicrobial Susceptibility Testing; Eleventh Informational Supplement. M100-S11, NCCLS, Wayne, PA, 2001.
- 6- Ahrabi, S. S., S. Kocagöz., T. Kocagöz, S. Ünal, A. Günalp. Evaluation of a New Rapid Antibacterial Susceptibility Test by comparison with disk diffusion method. American Society for Microbiology, 97th General Meeting. Miami Beach, Florida, U. S. A. May 4-8, 1997. C-487.

Manufacturer:

SALUBRIS, Inc.

Baldwin Park I,

12 Alfred Street, Suite 300

Woburn, MA 01801, U.S.A.

Tel: (800) 326-1484, Fax: (617) 249-0803

EUROPE Tel: +33 147 41 28 88

Website: www.salubrisinc.com

Authorised representative:

SALUBRIS

16, avenue de Beauval

92380 Garches - FRANCE

Tel: +33 1 47 41 28 88

Catalogue number:

MP 0250; MP 0250-12; MP 0250-15



QUICOLOR DISK DIFFUSION

INHIBITION ZONE INTERPRETIVE STANDARDS

Measured inhibition zone borders for resistance and susceptibility (x-y) in mm

Measured value lower than or equal to x means: Resistant;

Measured value higher than x and lower than y means : Intermediate;

Measured value higher than or equal to y means: Sensitive

<i>Antibacterial Agent</i>	<i>Disk content (µg)</i>	<i>Staphylococci</i>		<i>Enterobacteriaceae</i>		<i>Pseudomonadaceae (and other non-fermentative Gram (-) bacteria)</i>	
		<i>Mueller Hinton x-y (mm)</i>	<i>QuiColor x-y (mm)</i>	<i>Mueller Hinton x-y (mm)</i>	<i>QuiColor x-y (mm)</i>	<i>Mueller Hinton x-y (mm)</i>	<i>QuiColor x-y (mm)</i>
Amikacin	30	14-17	16-19	14-17	14-17	14-17	10-12
Amoxicillin/Clavulanate	20/10	19-20	19-20	13-18	13-17		
Ampicillin	10	28-29	22-23	13-17	14-16		
Ampicillin/Sulbactam	10/10	11-15	12-15	11-15	12-14		
Azlocillin	75					17-18	16-17
Aztreonam	30			15-22	12-17	15-22	12-17
Carbenicillin	100			19-23	15-20	13-17	13-15
Cefaclor	30	14-18	14-18	14-18	17-19		
Cefazolin	30	14-18	13-17	14-18	16-18		
Cefepime	30	14-18	15-17	14-18	16-18	14-18	14-16
Cefoperazone/Sulbactam	75/30			15-21	15-18	15-21	15-18
Ceftazidime	30	14-18	12-15	14-18	15-17	14-18	13-15
Ceftriaxone	30	13-21	13-17	13-21	15-21	13-21	12-14
Cefuroxime	30	14-23	14-18	14-23	15-20		
Ciprofloxacin	5	15-21	13-16	15-21	15-19	15-21	13-16
Clindamycin	2	14-21	14-20				
Erythromycin	15	13-23	17-20				
Fusidic acid	10	15-22	13-15				
Gentamicin	10	12-15	11-13	12-15	12-14	12-15	9-10
Imipenem	10			13-16	13-15	13-16	20-23
Meropenem	10			13-16	13-15	13-16	17-20
Netilmicin	30	12-15	11-13	12-15	13-15	12-15	9-11
Norfloxacin	10	12-17	12-15	12-17	13-15	12-17	11-13
Ofloxacin	5	12-16	14-16	12-16	13-16	12-16	11-14
Oxacillin	1	10-13	12-13				
Penicillin	10 units	28-29	23-24				
Piperacillin	100			17-21	15-17	17-18	15-16
Rifampin		16-20	16-20				
Teicoplanin	30	10-14	10-12				
Tetracycline	30	14-19	17-19	14-19	14-16		
Trimethoprim/Sulfamethoxazole	1.25/23.75	10-16	13-14	10-16	9-12		
Tobramycin	10	12-15	13-15	12-15	13-15	12-15	10-12
Vancomycin	30	9-12	9-12				